



# FORM DESIGNEHR RELEASE NOTES

VERSION 6.0.30

REQUIRES ECHOVANTAGE VERSION 3.24 OR HIGHER


## IMPROVEMENTS

### FDNEXT FORM PROPERTY, EV-385


A new form property to support opening the FDNext Client Allergies form from existing Form DesignEHR forms was added. When configured, the new Client Allergies form is opened in History mode. An existing allergy record can be selected, or use the blue plus “+” to add a new record. Existing allergy records created with the previous version of the *Allergies* form are opened in preview mode and have an upgrade option. See the EV v3.24 Release Notes for more information.

#### Configuring the FDNext Form Link

1. Open the form that requires the link to the New Allergies form in Form DesignEHR version 6.0.30.
2. Add a button.



3. With the newly added button selected (outlined in red), click into the Form Properties and make the following updates.




- a. Add a property *Name*, and *Hint Text* if desired.
- b. Add a *Caption*—the *Caption* text displays on the button/form link.
- c. Select **Form Link** in the *Button Type* dropdown listing.
- d. Select the *Form Link To Column* and *Form Link Value Column*. In this example,
  - i. The FHH form is linked to the *Client* column in the *dbo.ClientAllergies* table.
  - ii. The *Form Link Value Column* is the *Client* column in the *dbo.FamilyHealthHistories* table.
  - iii. This allows records for only the Client in focus to display when the form link is used.
- e. Ensure that the *Visible* checkbox is selected.

**\*\*Remember to save your updates.\*\***



### Using the New FDNNext Form Link

Once the form link is configured, open the form (FHH in this example) from the VHR Timeline.



When the newly created form link is selected, the New Allergies form opens in History mode. The history for references the FDNNext form's name, Client Allergies, not the caption used for the button.

