



# FORM DESIGNEHR RELEASE NOTES

VERSION 6.0.35


REQUIRES ECHOVANTAGE VERSION 3.32 OR HIGHER

## IMPROVEMENTS

### MASKED FIELDS, *EV-2248*

An issue was addressed where Form DesignEHR saved a masked field's mask characters in error when a user changed an existing masked field and then exited without saving. The masked values (\*\*\*\*) were saved in the database.

Form DesignEHR was updated to prevent this behavior by clearing the mask field when the user tabs or clicks into it and then setting the field to the same values if the user exits without making any changes.

Unselected Masked Field	
Selected Masked Field	